



1. Move all .nd2 into a new folder – name it whatever you like!
2. Run **batch-nd2-to-maxZ-sumZ-avgZ-plus-bmp.ijm** on the folder
 - a. Output: splits .nd2 into separate brightfield channel and three Z projections
 - b. Also an empty folder called “yeastspotter” that you need to populate by hand with segmented images
 - c. You can check the images by opening the Output-RGBbmp folder. Example below. From left to right: Average Z, Sum Z, Max Z of SEN555



3. Go to <http://yeastspotter.csb.utoronto.ca/> and upload 10 of the brightfield tifs at a time
 - a. I recommend you open a bunch of tabs & upload 10 brightfield images at a time into each new tab
4. Download the zip files when it is finished and extract them into the yeastspotter folder
 - a. Output: segmentation of your yeast based on the brightfield image. File names will correspond with the brightfield image name.
 - b. The only files in the yeastspotter folder should be the output images. Do not include the zip files or brightfield images in this folder.



5. Next, open and run the **quantify_perc_puncta_ver3_yeastspotter_detection.ijm** script in imageJ. You will be prompted to select the directory of Z projection you want to analyze

- a. I recommend maxZ for puncta; averageZ for fluorescence quantification if Z stacks are all different sizes (e.g average if you have 10 slices vs 30 slices in different samples), or sumZ if you want total cell fluorescence and all your Z stacks are the same size.
 - b. You will need to set the puncta threshold – this is the minimum value of intensity that is counted as a “puncta”. I recommend setting at 550 – this is what we used for the manual analysis before – but you might want to spot-check all your mutants and images to test this. Make sure you use the same threshold for ALL the genotypes/experiments you want to plot together! We should also report this value in the methods.
6. Compile results using the same R script as before!